


STANDARD OPERATING PROCEDURE			
	DNA Extraction from PBMCs (Qiagen Kit)		SOP #
			PL-SOP-003
	Originated by:	Ildiko Toth	Date:
	Laboratory:	Processing Laboratory	Pages:
	Approved by:	Alicja Trocha	25 April 17 1 of 2

I. PURPOSE:

The purpose of this procedure is to extract DNA from frozen PBMC samples.

II. SAFETY:

This protocol needs to be carried out in a biosafety cabinet in BSL 2+ Lab. Personnel are required to wear safety glasses, gloves and a disposable gown before initiating any work and a secondary pair of gloves over the first before entering the biosafety cabinet. Secondary pair of gloves **MUST** be removed every time while exiting from the hood even if it is for a short time.

III. REAGENTS:


Reagent	Vendor	Catalogue #
Dulbecco's Phosphate Buffered Saline	Sigma	D8537
QIAGEN protease*	Qiagen	19157
Buffer AL*	Qiagen	19075
QIAamp Spin Column*	Qiagen	cannot be ordered separately from kit
100% ethanol	Pharmco-AAPER	111ACS200
Buffer AW1*	Qiagen	19081
Buffer AW2*	Qiagen	19072
Buffer AE*	Qiagen	19077

*Reagents also part of "QIAamp DNA Blood Mini Kit", Cat. # 51104 (50 samples), 51106 (250 samples)

IV. PROCEDURE:

1. Thaw sample. Using approximately 10M cells for the extraction should produce a good yield. Spin down and resuspend cells in **200ul PBS**.
2. Add **200ul PBS/PBMC** mixture to **20ul QIAGEN protease** aliquot (found in -20 in Room 5225).
3. Add **200ul buffer AL** to the tube.
4. Vortex for 15 sec and incubate at 56° for 10 minutes. Centrifuge to remove condensation from the lid.
5. Add **200ul ethanol** (96-100%). Vortex and centrifuge.
6. Add mixture (about 650ul) to QIAamp Spin Column, in 2ml collection tube. Close cap and centrifuge.
7. Place column in clean collection tube and add **500ul buffer AW1**. Close cap and centrifuge.
8. Place column in clean collection tube and add **500ul buffer AW2**. Close cap and centrifuge at full speed (20,000 x g/14,000 rpm) for 3 minutes.
9. Place column in a new 1.5ml microcentrifuge tube. Add **100ul buffer AE** and incubate at room temperature overnight.

STANDARD OPERATING PROCEDURE

 Ragon Institute of MGH, MIT and Harvard	DNA Extraction from PBMCs (Qiagen Kit)		SOP #	PL-SOP-003
	Originated by:	Ildiko Toth	Date:	25 April 17
	Laboratory:	Processing Laboratory	Pages:	2 of 2
	Approved by:	Alicja Trocha		

10. The next day: centrifuge incubation sample.

11. Measure DNA concentration and aliquot samples for sequencing using nanodrop:

- 50ng/ul in 100 ul (for HLA typing)
- stock aliquot in 100 ul

V. REFERENCES/ ADDITIONS/ NOTES:

1. Set up 56° water bath prior to starting.
2. All centrifugation is done at 6,000 x g (8,000 rpm) for 1 min unless noted.