


STANDARD OPERATING PROCEDURE				
	<b>DNA Measurement using NanoDrop ND-1000 instrument</b>		<b>SOP #</b>	PL-SOP-006
	Originated by:	Ildiko Toth	Date:	27 April 17
	Laboratory:	Processing Laboratory	Pages:	1 of 2
	Approved by:	Alicja Trocha		

## I. PURPOSE:

The purpose of this procedure is DNA quantification using NanoDrop ND-1000 instrument.


## II. SAFETY:

This protocol needs to be carried out in a BSL1 Lab. Personnel are required to wear gloves and a lab coat before initiating any work.

## III. PROCEDURE:

### Wear Gloves.

1. Open program on computer.
2. Initialize spectrophotometer:
  - a. Load 2ul distilled water on measurement pedestal.
  - b. Lower sampling arm gently.
  - c. Click "OK".
3. Measure Blank:
  - a. Load 2ul blank buffer (the buffer, solvent, or carrier used in your sample).
  - b. Lower sampling arm gently.
  - c. Type BLANK into the Sample ID window.
  - d. Click "BLANK" (F3).
4. Ensure Quality of Blank:
  - a. Load 2ul of same blank buffer.
  - b. Lower sampling arm.
  - c. Click "MEASURE" (F1).
  - d. Check that the value in the concentration window is within -0.9 and +0.9.
  - e. If the value exceeds these limits, repeat the blanking process. If not, proceed to measure your DNA.
  - f. Repeat blanking process after every 10 samples.
5. Briefly vortex or mix DNA samples before loading on the NanoDrop.
6. Measure DNA concentration and purity:
  - a. Type in the Sample ID.
  - b. Load 2ul of DNA sample.
  - c. Lower sampling arm.
  - d. Click "MEASURE" (F1).
    - 260/280 window: a value of ~1.8 indicates pure DNA. Lower values indicate contaminants such as protein.
    - 260/230 window: a value of 1.8-2.2 indicates pure DNA. Lower values indicate contaminants

STANDARD OPERATING PROCEDURE				
	<b>DNA Measurement using NanoDrop ND-1000 instrument</b>		<b>SOP #</b>	PL-SOP-006
	Originated by:	Ildiko Toth	Date:	27 April 17
	Laboratory:	Processing Laboratory	Pages:	2 of 2
	Approved by:	Alicja Trocha		

- ng/ul window: displays DNA concentration.
7. Between each measurement, clean upper and lower measurement pedestals with de-ionized (miliQ) water on a kim-wipe.
  8. When you are finished measuring your samples clean the pedestals again and leave a dry, folded kim-wipe between the sampling arm and the lower measurement pedestal.
  9. Your measurements are automatically saved and can be found in: My Computer → Local Disk C → NanoDrop Data → Default → Nucleic Acids
    - a. Look for the date on which the measurements were taken.
  10. Remember to clean the area after your work is done

**IV. REFERENCES/ ADDITIONS/ NOTES:**

N/A