


| STANDARD OPERATING PROCEDURE  |                                 |                   |        |                      |
|---|---------------------------------|-------------------|--------|----------------------|
|  | T-cell Re-stimulation Procedure |                   | SOP#   | TC-SOP-006           |
|   | Originated by:                  | Alicja Trocha     | Date:  | 23<br>November<br>17 |
|   | Laboratory:                     | Walker Laboratory | Pages: | 1 of 2               |
|   | Approved by:                    | Alicja Trocha     |        |                      |

## I. PURPOSE

The purpose of this procedure is to outline the T-cell re-stimulation procedure performed in Ragon laboratories.

## II. SAFETY

This procedure needs to be carried out in the BSL2+ laboratory following all BL2+ regulations.

## III. SCOPE

This procedure applies to all employees, students, contractors and visitors that work on the CTL Feeding protocol.

## IV. REQUIREMENTS


Training to be obtained by qualified trainer or Subject Matter Expertise (SME).

## V. RESPONSIBILITIES

- A. The Ragon Institute **qualified trainers** are responsible for the overall implementation of this procedure and ensuring compliance and for periodic review of this procedure. Updates if any may be initiated by the qualified trainers or Subject Matter Expertise (SME).
- B. **All employees, students, contractors and visitors** are required to follow this procedure.

## VI. REAGENTS

| Item                                     | Manufacturer            | Order Number                  |
|--|-------------------------|-------------------------------|
| RPMI-1640                                | Sigma                   | R0883                         |
| Pen/Strep (5000 IU Pen/ 5000ug/mL Strep) | Mediatech               | 30-001-C1                     |
| L-glutamine (200mM; 29.2 mg/mL)          | Mediatech               | 25-002-C1                     |
| HEPES (1M; 238.3mg/mL)                   | Mediatech               | 25-060-C1                     |
| FBS, Heat-inactivated                    | Sigma                   | F4135 (a lot tested for here) |
| IL2                                      | Hoffmann-La Roche       | Ro 23-6019                    |
| PBS                                      | Sigma                   | D8537                         |
| PHA                                      | Fisher                  | R30852801                     |
| Anti-CD3 Antibody                        | NIH AIDS Reagent Rental |                               |
| 15 ml propylene tubes                    | Falcon                  | 352063                        |
| T25 flasks with 0.2 ul filter            | Corning                 | 3056                          |

| STANDARD OPERATING PROCEDURE  |                                 |                   |        |                      |
|---|---------------------------------|-------------------|--------|----------------------|
|  | T-cell Re-stimulation Procedure |                   | SOP#   | TC-SOP-006           |
|   | Originated by:                  | Alicja Trocha     | Date:  | 23<br>November<br>17 |
|   | Laboratory:                     | Walker Laboratory | Pages: | 2 of 2               |
|   | Approved by:                    | Alicja Trocha     |        |                      |

## VII. PROCEDURE

If you want to maintain the culture of T cells continuously for longer period, it is a good idea to restimulate every week by maintaining more than one batch of cells for each cell line. You will always keep older badges up so you can take 1-2M from those flasks for restimulation. Since clones grow differently restimulation to restimulation (due to different badges of fresh feeders), this will increase your chances of having enough cells for an experiment or for freezing.

Use Buffy coat cells (from Blood Bank at Main Campus) to restimulate or PBMC from HIV negative donor if smaller restimulation is planned (that can be requested through RAGON CLINICAL Platform. Follow protocol 00-03) for cloning use same day fresh blood/buffy coat.

**Note:** Always make fresh R10/50 for restimulation! Use fresh 12F6 ( anti CD3 ab) or PHA. (These should not be stored in the refrigerator for more than 2 weeks. Stock is kept in -80<sup>0</sup> freezer.)

1. Ficoll buffy coat or peripheral blood from HIV negative donor. See protocol .
2. After final wash, count cells and resuspend them at 20 million cells/ml of R10.
3. Irradiate at 3000 cGY in the XRAD irradiator on 9<sup>th</sup> floor room 917
4. Place 1-2 million CTLs in each T25 flask. Do not exceed 2 million CTLs per flask.
5. Add 20 million irradiated feeders in 20 ml of volume of R10/50 medium.
6. Add either 12F6 at 0.1  $\mu$ g/ml or PHA at 1:1000 dilutions.
7. Warning: do not restimulate more than 3-5 flasks at once! You can make 3-5 place them in the incubator and start new badge. This will ensure that the clones spend a minimal amount of time outside their desired incubator conditions. If you are a beginner, it's a good idea to do even fewer simultaneously because you are slower.
8. In keeping with that warning, make sure that your CTLs are not in the incubator which is opened frequently. This will help keep the incubator temperature and CO<sub>2</sub> levels as optimal as possible.
9. After CTLs have been established, they are fed twice a week (Monday /Friday) by medium exchange. Transfer flasks gently into the hood, being careful not to disturb the cells on the bottom. Aspirate about 7 ml of media off the top, and then add about 7 ml of fresh R10/50 to each flask. You can feed with media which is not older than 2 weeks.
10. Flasks are fed in this manner for several weeks or as long as the clone is growing.
11. Small numbers of cells can be removed and restimulated to make a fresh batch. Old cells will remain active as long as they are fed.
12. The recommended wait time for doing a Cr51 assay on restimulated CTLs is 7 days after restimulation. After that, the clone can be used as long as it is in good condition, which varies from clone to clone, batch to batch.
13. If testing specificity by Flow /Tetramer 12 days or more waiting period is recommended, unless you will use dead cell dye to sort out dead feeders.